

L11 ANSWER 1 OF 4 CA COPYRIGHT 2003 ACS
 AN 137:334779 CA
 TI FAST CARS: Engineering a laser spectroscopic technique for rapid
 identification of bacterial **spores**
 AU Scully, M. O.; Kattawar, G. W.; Lucht, R. P.; Opatrny, T.; Pilloff, H.;
 Rebane, A.; Sokolov, A. V.; Zubairy, M. S.
 CS Institute for Quantum Studies, Texas A and M University, College Station,
 TX, 77843, USA
 SO Proceedings of the National Academy of Sciences of the United States of
 America (2002), 99(17), 10994-11001
 CODEN: PNASA6; ISSN: 0027-8424
 PB National Academy of Sciences
 DT Journal
 LA English
 AB Airborne contaminants, e.g., bacterial **spores**, are usually
 analyzed by time-consuming microscopic, chem., and biol. assays. Current
 research into real-time laser spectroscopic detectors of such contaminants
 is based on e.g., resonance **fluorescence**. The present approach
 derives from recent expts. in which atoms and mols. are prepd. by one (or
 more) coherent laser(s) and probed by another set of lasers. However,
 generating and using maximally coherent oscillation in macromols. having
 an enormous no. of degrees of freedom is challenging. In particular, the
 short dephasing times and rapid internal conversion rates are major
 obstacles. However, adiabatic fast passage techniques and the ability to
 generate combs of phase-coherent femtosecond pulses provide tools for the
 generation and utilization of maximal quantum coherence in large mols. and
 biopolymers. We call this technique FAST CARS (femtosecond adaptive
 spectroscopic techniques for coherent anti-Stokes Raman spectroscopy), and
 the present article proposes and analyses ways in which it could be used
 to rapidly identify preselected mols. in real time.
 CC 9-5 (Biochemical Methods)
 ST FAST CARS engineering laser spectroscopic technique bacterial
spore
 IT Bacillus cereus
 Bacillus megaterium
Bacteria (Eubacteria)
 CARS Raman spectroscopy
 Laser spectroscopy
 Resonance **fluorescence**
Spore
 (FAST CARS for rapid identification of bacterial **spores**)
 IT Biopolymers
 Macromolecular compounds
 RL: ANT (Analyte); **ANST (Analytical study)**
 (FAST CARS for rapid identification of bacterial **spores**)
 IT DNA
 RL: BSU (Biological study, unclassified); BIOL (Biological study)
 (FAST CARS for rapid identification of bacterial **spores**)
 IT Proteins
 RL: BSU (Biological study, unclassified); BIOL (Biological study)
 (FAST CARS for rapid identification of bacterial **spores**)
 IT PNA
 RL: BSU (Biological study, unclassified); BIOL (Biological study)
 (FAST CARS for rapid identification of bacterial **spores**)
 IT 6893-30-7, Calcium **dipicolinate**
 RL: ARU (Analytical role, unclassified); **ANST (Analytical study)**
 (FAST CARS for rapid identification of bacterial **spores**)
 RE.CNT 66 THERE ARE 66 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L11 ANSWER 2 OF 4 CA COPYRIGHT 2003 ACS
 AN 134:233815 CA
 TI Physical perturbation for fluorescent characterization of microorganism

particles

AU Bronk, Burt V.; Shoaibi, Azadeh; Nudelman, Raphael; Akinyemi, Agnes N.
 CS AFRL/ECBC at U.S. Army ECBC, A.P.G., MD, 21010-5424, USA
 SO Proceedings of SPIE-The International Society for Optical Engineering
 (2000), 4036 (Chemical and Biological Sensing), 169-180
 CODEN: PSISDG; ISSN: 0277-786X

PB SPIE-The International Society for Optical Engineering
 DT Journal
 LA English

AB The motivation for using response to phys. perturbation to classify
 microparticles came from our previous expts. with **Dipicolinic**
 Acid (DPA). DPA as a calcium complex is a major component of bacterial
spores, constituting more than 5% of their dry wt. It is not
 commonly found in other natural products and therefore its presence is
 indicative of the presence of bacterial **spores**. Previous
 schemes utilizing the presence of DPA to detect these **spores**
 have relied on **fluorescence** which occurs when lanthanide metals
 (e.g., terbium) are added to a soln. where the presence of DPA is to be
 detd. We have recently demonstrated that changes in the
fluorescence of DPA can be stimulated without the addn. of such
 reagents. Thus after exposure to UV light, a substantial increase of
fluorescence emitted by DPA solns. with a peak at 410 nm occurs
 for excitation light with wavelength less than approx. 305 nm.

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10

ST fluorometry microorganism optical classification **dipicolinate**

IT Bacillus megaterium
 Bacillus subtilis
 Escherichia coli
 Fluorometry
 Microorganism
 (phys. perturbation for fluorescent characterization of microorganism
 particles)

IT 499-83-2, **Dipicolinic** Acid
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (phys. perturbation for fluorescent characterization of microorganism
 particles)

RE.CNT 11 THERE ARE 11 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L11 ANSWER 3 OF 4 CA COPYRIGHT 2003 ACS

AN 132:75561 CA

TI **Dipicolinic** acid (DPA) assay revisited and appraised for
spore detection

AU Hindle, Alistair A.; Hall, Elizabeth A. H.
 CS Inst. Biotechnol., University of Cambridge, Cambridge, CB2 1QT, UK
 SO Analyst (Cambridge, United Kingdom) (1999), 124(11), 1599-1604
 CODEN: ANALAO; ISSN: 0003-2654

PB Royal Society of Chemistry
 DT Journal
 LA English

AB Delayed gate **fluorescence** detection of **dipicolinic**
 acid (DPA), a universal and specific component of bacterial **spores**
 , has been appraised for use in a rapid anal. method for the detection of
 low concns. of bacterial **spores**. DPA was assayed by
 fluorimetric detection of its chelates with lanthanide metals. The
 influence of the choice and concn. of lanthanide and buffer ions on the
fluorescence assay was studied as well as the effects of pH and
 temp. The optimal system quantified the **fluorescence** of terbium
 monodipicolinate in a soln. of 10 .mu.M terbium chloride buffered with 1 M
 sodium acetate, pH 5.6 and had a detection limit of 2 nM DPA. This assay
 allowed the first real-time monitoring of the germination of bacterial
spores by continuously quantifying exuded DPA. A detection limit

of 104 *Bacillus subtilis* **spores** ml-1 was reached, representing a substantial improvement over previous rapid tests.

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10

ST **dipicolinate** DPA detn *Bacillus* **spore** fluorometry

IT *Bacillus subtilis*
 Fluorometry
 Solvent effect
Spore
 (**dipicolinic** acid (DPA) assay for **spore** detection)

IT 499-83-2, **Dipicolinic** acid
 RL: ANT (Analyte); BOC (Biological occurrence); BSU (Biological study, unclassified); **ANST (Analytical study)**; BIOL (Biological study); OCCU (Occurrence)
 (**dipicolinic** acid (DPA) assay for **spore** detection)

IT 13759-92-7, Europium chloride hexahydrate 13798-24-8, Terbium chloride hexahydrate 15059-52-6, Dysprosium chloride hexahydrate
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES (Uses)
 (**dipicolinic** acid (DPA) assay for **spore** detection)

RE.CNT 51 THERE ARE 51 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L11 ANSWER 4 OF 4 CA COPYRIGHT 2003 ACS

AN 131:141627 CA

TI **Fluorescence** of **dipicolinic** acid as a possible component of the observed UV emission spectra of bacterial **spores**

AU Nudelman, Raphael; Feay, Nicole; Hirsch, Mathew; Efrima, Shlomo; Bronk, Burt

CS Mantech Environmental Technology Inc., USA

SO Proceedings of SPIE-The International Society for Optical Engineering (1999), 3533 (Air Monitoring and Detection of Chemical and Biological Agents), 190-195
 CODEN: PSISDG; ISSN: 0277-786X

PB SPIE-The International Society for Optical Engineering

DT Journal

LA English

AB **Dipicolinic** acid (DPA) and the Ca²⁺ complex of DPA (CaDPA) are well-known and are major chem. components of bacterial **spores**. DPA's native **fluorescence** is very weak and is thought to be completely masked by the **fluorescence** of tryptophan when this compd. is present. Thus **fluorescence** related to DPA in **spores** is assumed by many authors to be completely absent. We show that the **fluorescence** of CaDPA is substantial for excitation between about 290 nm and 310 nm with emission peaking near 400 nm. This emission is at the long wavelength tail for emission from tryptophan. We examine whether the emission of CaDPA could contribute to the total emission spectrum when bacterial **spores** are present in an aerosol, for excitation wavelengths in the neighborhood of 310 nm. In this report we present measurements of **fluorescence** excitation and emission for CaDPA and compare them with that of DPA and tryptophan.

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 4, 10

ST bacterium **spore** **dipicolinate** fluorometry

IT **Bacteria** (Eubacteria)
 Environmental analysis
 Fluorometry
Spore
 (**fluorescence** of **dipicolinic** acid as a possible component of obsd. UV emission spectra of bacterial **spores**)

IT 73-22-3, L-Tryptophan, analysis 499-83-2, **Dipicolinic** acid 499-83-2D, **Dipicolinic** acid, calcium complex 7440-70-2D, Calcium, complex with **dipicolinic** acid, analysis
 RL: ANT (Analyte); BOC (Biological occurrence); BSU (Biological study,

unclassified); **ANST** (**Analytical study**); **BIOL** (**Biological study**);
OCCU (**Occurrence**)

(**fluorescence** of **dipicolinic acid** as a possible
component of obsd. UV emission spectra of bacterial **spores**;

RE.CNT 15 THERE ARE 15 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

=>

L19 ANSWER 4 OF 4 CA COPYRIGHT 2003 ACS
 AN 123:328802 CA
 TI Native **fluorescence** detection and spectral differentiation of
 peptides containing tryptophan and tyrosine in capillary electrophoresis
 AU Timperman, Aaron T.; Oldenburg, Kurt E.; Sweedler, Jonathan V.
 CS Department of Chemistry, University of Illinois, Urbana, IL, 61801, USA
 SO Analytical Chemistry (1995), 67(19), 3421-6
 CODEN: ANCHAM; ISSN: 0003-2700
 PB American Chemical Society
 DT Journal
 LA English
 AB A native **fluorescence** detection system for capillary
 electrophoresis is described that achieves low attomole detection limits
 and simultaneous acquisition of complete **fluorescence**
emission spectra. The system is designed for detection of
 peptides through the **intrinsic fluorescence** of
 tryptophan and tyrosine residues. The detection system employs a
 frequency doubled krypton laser operating at 284 nm for **excitation**
 , a sheath flow cell, a reflective f/1.2 microscope objective, an imaging
 spectrograph, and a CCD detector. The detection capabilities were
 characterized with tryptophan and tyrosine, which have limits of detection
 (3.sigma.) of 2 .times. 10-10 and, 2 .times. 10-8 M, resp. Acquisition of
 the **fluorescence emission** spectrum provides the
 ability to distinguish three classes of peptides: those that contain
 tryptophan, those that contain tyrosine, and those that contain both
 tryptophan and tyrosine.
 CC 80-2 (Organic Analytical Chemistry)
 Section cross-reference(s): 9
 ST native **fluorescence** detection tryptophan tyrosine peptide;
 electrophoresis **fluorescence** detection tryptophan tyrosine
 peptide
 IT Fluorometers
 (for native **fluorescence** detection and spectral
 differentiation of peptides contg. tryptophan and tyrosine in capillary
 electrophoresis)
 IT Peptides, analysis
 RL: ANT (Analyte); **ANST (Analytical study)**
 (native **fluorescence** detection and spectral differentiation
 of peptides contg. tryptophan and tyrosine in capillary
 electrophoresis)
 IT Electrophoresis and Ionophoresis
 (capillary, detectors, for native **fluorescence** detection and
 spectral differentiation of peptides contg. tryptophan and tyrosine)
 IT **Microbial** hormones and pheromones
 RL: ANT (Analyte); **ANST (Analytical study)**
 (.alpha.-factor, native **fluorescence** detection and spectral
 differentiation of peptides contg. tryptophan and tyrosine in capillary
 electrophoresis)
 IT 60-18-4D, Tyrosine, peptides 73-22-3, Tryptophan, analysis 73-22-3D,
 Tryptophan, peptides 21778-69-8 58822-25-6, Leucine-enkephalin
 59401-28-4, .alpha.-Mating factor (yeast) 65418-88-4, .alpha.-1-Mating
 factor 89911-64-8, Cholecystokinin(26-31) 98395-75-6, Neuromedin U 8
 (pig spinal cord)
 RL: ANT (Analyte); **ANST (Analytical study)**
 (native **fluorescence** detection and spectral differentiation
 of peptides contg. tryptophan and tyrosine in capillary
 electrophoresis)

=>

L17 ANSWER 5 OF 179 CA COPYRIGHT 2003 ACS
 AN 138:85885 CA
 TI Comparative intrinsic and enhanced total photoluminescence of endospore material
 AU Anderson, John E.; Webb, Stanley Thomas; Fischer, Robert L.; Kester, Karen; Smith, Clint
 CS U.S. Army Topographic Engineering Center, Alexandria, VA, 22315, USA
 SO Proceedings of SPIE-The International Society for Optical Engineering (2002), 4576 (Advanced Environmental Sensing Technology II), 27-31
 CODEN: PSISDG; ISSN: 0277-786X
 PB SPIE-The International Society for Optical Engineering
 DT Journal
 LA English
 AB Two techniques are compared using total luminescence spectroscopy to detect endospore material in preps. equiv. to $3.0 \times 10^5/\text{mL}$ **spores**. The first method applied intrinsic, steady-state photoluminescence for detection. The second approach using a binding fluorochrome derived from 4-p-dimethylaminostyrylpyridinium (DASP) to signal the presence of **spore** material. Comparative **fluorescence emission** signatures (excited at 469 nm) showed greater calibrated signal recovery (4×10^6 cps) for **spore** material at longer wavelengths using DASP. The intrinsic **fluorescence emission** of endospores (excited at 346 nm) occurred at shorter wavelengths and showed a reduced calibrated intensity (1.4×10^5 counts per s, cps). One major advantage of DASP appears to be its longer wavelength **excitation** (469 nm) that is out of the range of assocd. biol. materials that compete for absorption at shorter UV wavelengths.
 CC 9-5 (Biochemical Methods)
 ST photoluminescence spectroscopy endospore
 IT Luminescence spectroscopy
 (comparative intrinsic and enhanced total photoluminescence of endospore material using dimethylaminostyrylpyridinium fluorochrome)
 IT Fluorescent dyes
 (dimethylaminostyrylpyridinium deriv.; comparative intrinsic and enhanced total photoluminescence of endospore material using dimethylaminostyrylpyridinium fluorochrome)
 IT **Spore**
 (endospore; comparative intrinsic and enhanced total photoluminescence of endospore material using dimethylaminostyrylpyridinium fluorochrome)
 IT 186659-60-9
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES (Uses)
 (fluorochrome derived from; comparative intrinsic and enhanced total photoluminescence of endospore material using dimethylaminostyrylpyridinium fluorochrome)
 RE.CNT 10 THERE ARE 10 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 12 OF 179 CA COPYRIGHT 2003 ACS
 AN 137:286140 CA
 TI Optical structure for multi-photon **excitation** and the use thereof
 IN Duveneck, Gert L.; Bopp, Martin A.; Pawlak, Michael; Ehrat, Markus; Marowsky, Gerd
 PA Zeptosens A.-G., Switz.
 SO PCT Int. Appl., 76 pp.
 CODEN: PIXXD2
 DT Patent
 LA German
 FAN.CNT 1

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
-----	-----	-----	-----	-----

PI WO 2002079765 A2 20021010 WO 2002-EP2958 20020318
 WO 2002079765 A3 20030130
 W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN,
 CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,
 GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR,
 LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT,
 RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US,
 UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM
 RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH,
 CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR,
 BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG

PRAI CH 2001-617 A 20010402
 CH 2001-689 A 20010412

AB Optical structures comprising an optical waveguide with a waveguiding layer that is transparent at .gtoreq.1 **excitation** wavelength are described in which the intensity of **excitation** light that is input into and guided via the waveguiding layer is sufficiently high to excite mols. or groups of mols. at a surface of the layer or within 200 nm of the layer by means of multi-photon **excitation**, preferably two-photon **excitation**. Optionally, an adhesive agent may be used to immobilize biochem. or biol. mols., or other materials, for use in detecting or detg. analytes in a sample. Optical systems for carrying out multiphoton **excitation**, and methods for anal. (e.g., luminescence **excitation** and to the luminescent detection of one or several analytes) using the systems are also described. The systems may also be used to form optical tweezers. Methods for luminescence and **fluorescence** anal. of biomols., including autofluorescence of nucleic acids, is emphasized.

IC ICM G01N021-77
 ICS G01N021-55; G01N021-64

CC 73-10 (Optical, Electron, and Mass Spectroscopy and Other Related Properties)
 Section cross-reference(s): 9, 17, 64, 74

ST waveguide multiphoton **excitation** system; luminescence analysis
 waveguide multiphoton **excitation** system; **fluorescence**
 analysis waveguide multiphoton **excitation** system

IT Prion proteins
 RL: ANT (Analyte); ANST (Analytical study)
 (PrPSc, assay for; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT **Bacteria** (Eubacteria)
 Pathogen
 Salmonella
 (assay for; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT **Fluorescence**
 (autofluorescence; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Disease, plant
 (diagnosis of; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Immunoassay
 (**fluorescence**; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Immunoassay
 (immunofluorometric; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Photoexcitation

(multiphoton; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Polyesters, uses
 RL: DEV (Device component use); USES (Uses)
 (thio-; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Photoexcitation
 (two-photon; waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Blood analysis
 Combinatorial chemistry
 Diagnosis
 Drug design
 Drug screening
 Egg yolk
Fluorescence
 Fluorometers
 Fluorometry
 Food analysis
 Gene expression profiles
 Luminescence
 Optical waveguides
 Pharmaceutical analysis
 Plant analysis
 Surface waters
 Urine analysis
 (waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Agglutinins and Lectins
 Carbohydrates, analysis
 DNA
 Enzymes, analysis
 Oligonucleotides
 PNA
 Receptors
 PL: ANT (Analyte); **ANST (Analytical study)**
 (waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Acrylic polymers, uses
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Polyamides, uses
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Polycarbonates, uses
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Polyesters, uses
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods of multiphoton **excitation** using them and use of the structures in luminescent anal.)

IT Polyimides, uses

RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods
 of multiphoton **excitation** using them and use of the
 structures in luminescent anal.)

IT Polythiophenylenes
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods
 of multiphoton **excitation** using them and use of the
 structures in luminescent anal.)

IT Polyurethanes, uses
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods
 of multiphoton **excitation** using them and use of the
 structures in luminescent anal.)

IT 71-00-1D, L-Histidine, oligomers
 RL: ANT (Analyte); **ANST (Analytical study)**
 (waveguide structures for multiphoton **excitation** and methods
 of multiphoton **excitation** using them and use of the
 structures in anal.)

IT 1313-96-8, Niobium oxide 1314-13-2, Zinc oxide, uses 1314-23-4,
 Zirconium oxide, uses 1314-61-0, Tantalum oxide 9002-88-4,
 Polyethylene 9003-01-4, Polyacrylic acid 9003-07-0, Polypropylene
 9003-53-6, Polystyrene 9011-14-7, Polymethylmethacrylate 12055-23-1,
 Hafnium oxide 13463-67-7, Titanium oxide, uses 25038-59-9,
 Polyethylene terephthalate, uses
 RL: DEV (Device component use); USES (Uses)
 (waveguide structures for multiphoton **excitation** and methods
 of multiphoton **excitation** using them and use of the
 structures in luminescent anal.)

L17 ANSWER 17 OF 179 CA COPYRIGHT 2003 ACS

AN 136:382505 CA

TI Device for monitoring cells

IN Pitner, J. Bruce; Hemperly, John Jacob; Guarino, Richard D.; Wodnicka,
 Magdalena; Stitt, David T.; Burrell, Gregory J.; Foley, Timothy G., Jr.;
 Beaty, Patrick Shawn

PA Becton, Dickinson and Company, USA

SO U.S., 42 pp., Cont.-in-part of U.S. Ser. No. 715,557.

CODEN: USXXAM

DT Patent

LA English

FAN.CNT 3

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	US 6395506	B1	20020528	US 1999-342720	19990629
	EP 509791	A1	19921021	EP 1992-303391	19920415
	EP 509791	B1	19960703		
	R: DE, FR, GB, IT				
	CA 2066329	AA	19921019	CA 1992-2066329	19920416
	JP 05137596	A2	19930601	JP 1992-98368	19920418
	JP 07073510	B4	19950809		
	US 2002192636	A1	20021219	US 2002-109475	20020328
	US 2002155424	A1	20021024	US 2002-116777	20020404
PRAI	US 1991-687359	B1	19910418		
	US 1993-25899	A2	19930303		
	US 1996-715557	A2	19960918		
	US 1999-342720	A2	19990629		
	US 2000-642504	A2	20000818		
	US 2001-966505	A2	20010928		

AB The present invention relates to methods for detection and evaluation of
 metabolic activity of eukaryotic and/or prokaryotic cells based upon their
 ability to consume dissolved oxygen. The methods utilize a luminescence
 detection system which makes use of the sensitivity of the luminescent
emission of certain compds. to the presence of oxygen, which

quenches (diminishes) the compd.'s luminescent **emission** in a concn. dependent manner. Respiring eukaryotic and/or prokaryotic cells will affect the oxygen concn. of a liq. medium in which they are immersed. Thus, this invention provides a convenient system to gather information on the presence, identification, quantification and cytotoxic activity of eukaryotic and/or prokaryotic cells by detg. their effect on the oxygen concn. of the media in which they are present.

IC ICM C12Q001-18
NCL 435032000
CC 9-1 (Biochemical Methods)
Section cross-reference(s): 1, 4
ST device monitoring cell
IT Plates
(Microtitration; device for monitoring cells)
IT Analytical apparatus
Antibiotics
Bacteria (Eubacteria)
Biological materials
Blood
Blood serum
Cell
Cell proliferation
Chemicals
Coating materials
Composition
Concentration (condition)
Culture media
Cytotoxicity
Drugs
Escherichia coli
Eukaryota
Extracellular matrix
Fluorescence quenching
Impermeability
Insecta
Light
Liquids
Luminescence
Luminescence quenching
Luminescence spectroscopy
Luminescent substances
Mathematical methods
Metabolism
Microorganism
Molecules
Particles
Permeability
Prokaryote
Pseudomonas aeruginosa
Radiation
Reducing agents
Respiration, animal
Respiration, **microbial**
Sensors
Solutes
Wavelength
Wetting
Yeast
(device for monitoring cells)
IT Toxins
RL: ADV (Adverse effect, including toxicity); BIOL (Biological study)
(device for monitoring cells)
IT Reagents
RL: ARG (Analytical reagent use); **ANST** (Analytical study); USES

(Uses)
(device for monitoring cells)

IT Plastics, analysis
PL: ARU (Analytical role, unclassified); **ANST (Analytical study)**
(device for monitoring cells)

IT Pubber, analysis
PL: ARU (Analytical role, unclassified); **ANST (Analytical study)**
(device for monitoring cells)

IT Silicone rubber, analysis
PL: ARU (Analytical role, unclassified); **ANST (Analytical study)**
(device for monitoring cells)

IT Growth factors, animal
PL: BSU (Biological study, unclassified); BIOL (Biological study)
(device for monitoring cells)

IT Collagens, biological studies
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)
(device for monitoring cells)

IT Entactin
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)
(device for monitoring cells)

IT Laminins
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)
(device for monitoring cells)

IT Proteoglycans, biological studies
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)
(heparitin sulfate-contg.; device for monitoring cells)

IT Optical detectors
(luminescence; device for monitoring cells)

IT Animal cell
(mammal; device for monitoring cells)

IT Amino acids, biological studies
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)
(nonessential; device for monitoring cells)

IT Collagens, biological studies
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)
(type IV; device for monitoring cells)

IT 1499-10-1, 9,10-Diphenylanthracene 15158-62-0D, Tris-2,2'-
bipyridylruthenium (II), salts 36309-88-3, Tris-4,7-diphenyl-1,10-
phenanthroline ruthenium (II) chloride 50525-27-4, Tris-2,2'-
bipyridylruthenium (II) chloride hexahydrate. 63373-04-6D,
Tris-4,7-diphenyl-1,10-phenanthroline ruthenium (II), salts
PL: APG (Analytical reagent use); **ANST (Analytical study)**; USES
(Uses)
(device for monitoring cells)

IT 7631-86-9, Silica, analysis
PL: ARU (Analytical role, unclassified); **ANST (Analytical study)**
(device for monitoring cells)

IT 59-05-2, Methotrexate 151-21-3, Sodium dodecyl sulfate, biological
studies 865-21-4, Vinblastine 7757-83-7, Sodium Sulfite 7782-44-7,
Oxygen, biological studies 26628-22-8, Sodium Azide 35607-66-0,
Cefoxitin 55268-75-2, Cefuroxime 85721-33-1, Ciprofloxacin
PL: BSU (Biological study, unclassified); BIOL (Biological study)
(device for monitoring cells)

IT 57-92-1, Streptomycin, biological studies 113-24-6, Sodium pyruvate
1397-89-3, Fungizone 1406-05-9, Penicillin 119978-18-6, Matrigel
141907-41-7, Matrix metalloproteinase
PL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)

(device for monitoring cells)

RE.CNT 9 THERE ARE 9 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 18 OF 179 CA COPYRIGHT 2003 ACS

AN 136:337350 CA

TI Method for detecting living cells in medium and measuring pH of medium

IN Kawasaki, Yukishige; Tsuji, Takashi; Kurane, Ryuichiro

PA Sangyo Gijutsu Sogo Kenkyusho, Japan; Bioindustry Association; Mitsubishi Chemical Corp.

SO Jpn. Kokai Tokkyo Koho, 12 pp.

CODEN: JKXXAF

DT Patent

LA Japanese

FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	JP 2002125696	A2	20020508	JP 2000-319247	20001019
PRAI	JP 2000-319247		20001019		

AB A convenient method is provided for simultaneously performing the detection of living cells (e.g, microorganism) in a medium and the measurement of the medium pH. The method comprises a step for adding a fluorescent enzyme substrate (e.g., 5-carboxyfluorescein diacetate acetoxymethyl ester, 5-(and 6-)carboxyfluorescein diacetate, 2',7'-bis-(2-carboxyethyl)-5-(and 6-)carboxyfluorescein acetoxymethyl ester, 5-sulfofluorescein diacetate) to the medium contg. living cells, a step for irradiating two kinds of **excitation** light with different wavelength to the medium, a step for detecting the living cells from the **fluorescence** intensity obtained, and a step for measuring the pH of the medium by calcg. the intensity ratio between the **fluorescence** generated by two kinds of **excitation** light.

IC ICM C12Q001-06

ICS C12M001-34; C12Q001-34; G01N021-64; G01N021-80; G01N033-84

CC 9-5 (Biochemical Methods)

Section cross-reference(s): 10

ST living cell microorganism medium pH fluorometry

IT Optical filters

(bandpass; method for detecting living cells in medium and measuring pH of medium)

IT Fluorescent substances

(enzyme substrate; method for detecting living cells in medium and measuring pH of medium)

IT Light

(**excitation**; method for detecting living cells in medium and measuring pH of medium)

IT Cell

(living; method for detecting living cells in medium and measuring pH of medium)

IT Aeolosoma

Bacteria (Eubacteria)

Chlorella

Culture media

Cyclidium (protozoan)

Fluorescence microscopy

Fluorometry

Lepadella

Microorganism

Philodina

Scenedesmus

Schizothrix

Test kits

Wavelength

pH

(method for detecting living cells in medium and measuring pH of

medium)

IT 116723-31-0 118493-69-9 124387-19-5 124412-00-6
 RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
 (method for detecting living cells in medium and measuring pH of
 medium)

L17 ANSWER 20 OF 179 CA COPYRIGHT 2003 ACS

AN 136:337178 CA

TI Spectrofluorometric system devised so as to permit the synchronous
 measurement of both the bacterial bioluminescence and its related
fluorescence emission

AU Karatani, Hajime; Furuta, Kenji; Hirayama, Satoshi

CS Dept. of Polymer Science and Engineering, Kyoto Institute of Technology,
 Japan

SO Proceedings of SPIE-The International Society for Optical Engineering
 (2001), 4252 (Advances in Fluorescence Sensing Technology V), 88-96
 CODEN: PSISDG; ISSN: 0277-786X

PB SPIE-The International Society for Optical Engineering

DT Journal

LA English

AB A novel spectrofluorometer equipped with two photo-detectors has been
 developed based on a conventional spectrofluorometer with a view to
 measuring both bioluminescence and its related **fluorescence**
 originating from bacterial luciferase reaction. The configuration of the
 two monochromator-photomultiplier systems is to be opposite to each other.
 First the capability of the created system was evaluated using the
 peroxidase-catalyzed luminol chemiluminescence. The light
emission signals as a function of time evidently showed that
 intensities of chemiluminescence generated by the reaction and of
fluorescence elicited from the photoexcited reaction product vary
 inversely. This reflects the feature of the luminol reaction. Any
 interference from either **excitation** light or
fluorescence emission is absent in the detection of
 chem. initiated light. Subsequently, the evaluation of the
 spectrofluorometer was made on the luciferase reactions under the various
 conditions. The observation of signals for bioluminescence and
fluorescence from the luciferase reactions under the various
 conditions. The observation of signals for bioluminescence and
fluorescence from the luciferase reaction was established to be
 useful in studying the time-dependent behaviors of the fluorescent
 substrates and/or products as well as a primary emitter formed in the
 luciferase reaction. From these evaluations, the developed
 spectrofluorometer has proved to be profitable to study bacterial
 bioluminescence.

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10

ST spectrofluorometry measurement **bacteria** bioluminescence
fluorescence

IT **Bacteria** (Eubacteria)
Fluorescence
 Fluorometry
 Luminescence, bioluminescence
 (spectrofluorometric system for synchronous measurement of both
 bacterial bioluminescence and related **fluorescence**
emission)

IT 9014-00-0, Luciferase 39346-42-4, FMN reductase
 RL: ANT (Analyte); **ANST (Analytical study)**
 (spectrofluorometric system for synchronous measurement of both
 bacterial bioluminescence and related **fluorescence**
emission)

RE.CNT 12 THERE ARE 12 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 33 OF 179 CA COPYRIGHT 2003 ACS
 AN 135:192324 CA
 TI Ultraviolet **fluorescence** imaging applications
 AU Hill, Ralph H., Jr.; Angell, Peter
 CS Instrumentation and Space Research Division, Southwest Research Institute,
 San Antonio, TX, USA
 SO AT-PROCESS (2000), 5(3,4), 108-114
 CODEN: APJCFR; ISSN: 1077-419X
 PB InfoScience Services
 DT Journal
 LA English
 AB The UV **fluorescence** of arom. amino acids in **microbial**
 biofilms can be used to det. the biomass formed in corrosion pits on metal
 surfaces. This information is important in establishing the relationship
 between **bacteria** and corrosion; i.e., which comes first, the
 corrosion pits or the biomass. One specific amino acid that has been used
 in past studies to indicate biomass is tryptophane. Tryptophane
 fluoresces in the near UV region; this **fluorescence** can be used
 to quantify the amt. of tryptophane and biomass, present on the metal
 surfaces. Under Southwest Research Institute (SwRI) internal research
 funding, a radiometrically calibrated UV-imaging system has been
 previously developed. This system was developed to image tactical missile
 plumes in the solar-blind UV region of the spectrum near 260 nm. The
 basic building block is a microchannel plate-intensified charge-coupled
 device (MCP-CCD) camera. For the project presented in this paper, a 325
 nm helium-cadmium laser was used as an **excitation** source and the
 camera filtering converted to the near-UV wavelength region. Research was
 conducted to measure the **fluorescence** from the arom. amino acid
 tryptophane under various dilns. in water. Fluorescent images were also
 recorded from the Oceanospirillum **bacteria** on a copper coupon;
 this **bacteria** was originally isolated from copper on Navy
 platforms by Naval Research Lab. personnel. Other applications will also
 be mentioned.
 CC 9-5 (Biochemical Methods)
 ST UV **fluorescence** imaging
 IT Lasers
 (Helium-cadmium; UV **fluorescence** imaging applications)
 IT **Bacteria** (Eubacteria)
 Biofilms (**microbial** reactors)
 Biomass
 CCD cameras
 Corrosion
 Dilution
Fluorescence
 Interface
 Oceanospirillum
 Plates
 (UV **fluorescence** imaging applications)
 IT Metals, biological studies
 RL: BSU (Biological study, unclassified); BIOL (Biological study)
 (UV **fluorescence** imaging applications)
 IT Imaging
 (UV **fluorescence**; UV **fluorescence** imaging
 applications)
 IT Amino acids, analysis
 RL: ANT (Analyte); BSU (Biological study, unclassified); ANST
 (**Analytical study**); BIOL (Biological study)
 (arom.; UV **fluorescence** imaging applications)
 IT Wavelength
 (near-UV; UV **fluorescence** imaging applications)
 IT 73-22-3, Tryptophane, analysis
 RL: ANT (Analyte); BSU (Biological study, unclassified); ANST
 (**Analytical study**); BIOL (Biological study)
 (UV **fluorescence** imaging applications)

IT 7440-50-8, Copper, biological studies 7732-18-5, Water, biological studies
RL: BSU (Biological study, unclassified); BIOL (Biological study);
(UV **fluorescence** imaging applications)

RE.CNT 9 THERE ARE 9 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 36 OF 179 CA COPYRIGHT 2003 ACS

AN 135:18679 CA

TI Association of bright greenish yellow **fluorescence** with aflatoxin production in cereals

AU Ijaz, Nusrat; Salah-ud-Din; Yasin, M.

CS Biotechnology and Food Research Centre, Pakistan Council of Scientific and Industrial Research Laboratories, Lahore, Pak.

SO Pakistan Journal of Science (2000), 52(1-2), 47-52
CODEN: PAJSAS; ISSN: 0030-9877

PB Pakistan Association for the Advancement of Science

DT Journal

LA English

AB Bright greenish yellow **fluorescence** (BGYF) was studied in cereals, after inoculating those with **spores** of toxic strain of *Aspergillus flavus*. Pos. relationship between mold growth, kojic acid prodn., BGYF **emission** and aflatoxin contents was established. BGYF units increased with an increase in aflatoxin prodn. Profused mold growth resulted in improved synthesis of Kojic acid, the chem. which is necessary to produce BGYF.

CC 17-1 (Food and Feed Chemistry)

ST cereal aflatoxin **fluorescence** fluorimeter

IT *Aspergillus flavus*

Cereal (grain)

Fluorescence

Fluorometers

Rice (*Oryza sativa*)

Sorghum

Wheat

(assocn. of bright greenish yellow **fluorescence** with aflatoxin prodn. in cereals)

IT Aflatoxins

RL: ANT (Analyte); BOC (Biological occurrence); BSU (Biological study, unclassified); **ANST (Analytical study)**; BIOL (Biological study); OCCU (Occurrence)

(assocn. of bright greenish yellow **fluorescence** with aflatoxin prodn. in cereals)

IT 501-30-4, Kojic acid 9003-99-0, Peroxidase

RL: BOC (Biological occurrence); BSU (Biological study, unclassified); BIOL (Biological study); OCCU (Occurrence)

(assocn. of bright greenish yellow **fluorescence** with aflatoxin prodn. in cereals)

RE.CNT 17 THERE ARE 17 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 42 OF 179 CA COPYRIGHT 2003 ACS

AN 134:233815 CA

TI Physical perturbation for fluorescent characterization of microorganism particles

AU Bronk, Burt V.; Shoaibi, Azadeh; Nudelman, Raphael; Akinyemi, Agnes N.

CS AFRL/ECBC at U.S. Army ECBC, A.P.G., MD, 21010-5424, USA

SO Proceedings of SPIE-The International Society for Optical Engineering (2000), 4036 (Chemical and Biological Sensing), 169-180

CODEN: PSISDG; ISSN: 0277-786X

PB SPIE-The International Society for Optical Engineering

DT Journal

LA English

AB The motivation for using response to phys. perturbation to classify

microparticles came from our previous expts. with Dipicolinic Acid (DPA). DPA as a calcium complex is a major component of bacterial **spores**, constituting more than 5% of their dry wt. It is not commonly found in other natural products and therefore its presence is indicative of the presence of bacterial **spores**. Previous schemes utilizing the presence of DPA to detect these **spores** have relied on **fluorescence** which occurs when lanthanide metals (e.g., terbium) are added to a soln. where the presence of DPA is to be detd. We have recently demonstrated that changes in the **fluorescence** of DPA can be stimulated without the addn. of such reagents. Thus after exposure to UV light, a substantial increase of **fluorescence** emitted by DPA solns. with a peak at 410 nm occurs for **excitation** light with wavelength less than approx. 305 nm.

CC 9-5 (Biochemical Methods)
Section cross-reference(s): 10
ST fluorometry microorganism optical classification dipicolinate
IT Bacillus megaterium
Bacillus subtilis
Escherichia coli
Fluorometry
Microorganism
(phys. perturbation for fluorescent characterization of microorganism particles)
IT 499-83-2, Dipicolinic Acid
RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
(Uses)
(phys. perturbation for fluorescent characterization of microorganism particles)
RE.CNT 11 THERE ARE 11 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 46 OF 179 CA COPYRIGHT 2003 ACS
AN 134:128009 CA
TI Sensitivity of detection of **bacteria** with fluorescent and luminescent phenotypes using different instruments
AU Brovko, Liubov Yu.; Griffiths, Mansel W.
CS Food Sci. Dep., Univ. of Guelph, Guelph, ON, Can.
SO Proceedings of SPIE-The International Society for Optical Engineering (2000), 3921(Optical Diagnostics of Living Cells III), 147-156
CODEN: PSISDG; ISSN: 0277-786X
PB SPIE-The International Society for Optical Engineering
DT Journal
LA English
AB The problem of bacterial enumeration in different samples is of great importance in many fields of research. Construction of recombinant fluorescent and luminescent **bacteria** that can be easily detected by nondestructive instrumental methods provides us with an opportunity to monitor **bacteria** in a wide variety of clin., environmental and food samples in real time. Three different labels were employed: Green Fluorescent Protein (GFP), Bacterial luciferase (BL) and Firefly Luciferase (FFL). Both plasmid and chromosomal transformants of different strains of E. coli, P. putida and S. enteritidis were used. For the detection of the in vivo GFP the Shimadzu RF 540 spectrofluorimeter, Labsystems FL- 500 plate fluorometer and Night Owl LB 98 CCD-camera from EG and G Berthold supplied with **excitation** light source and proper spectral filters both in macroscopic and microscopic mode were used. For the detection of in vivo luminescence of BL and FFL, tube luminometer BG-P from GEM Biomedical Inc., luminometric plate reader from BioOrbit, BIQ Bioview CCD camera from Cambridge Imaging Ltd and Night Owl LB 98 CCD camera both in macroscopic and microscopic mode were used. The expression levels of the labels, their stability, stability of the signal and detection limits of tagged **bacteria** were investigated. The detection limits for GFP tagged **bacteria** were 5×10^4 - 6×10^6 , for BL tagged **bacteria** 5×10^2 - 2×10^5 , and for FFL tagged

bacteria - 4×10^3 - 10^6 CFU/mL, depending on the instrument used.
Single **bacteria** could be detected with the help of the Night Owl
in the microscopic mode.

CC 9-5 (Biochemical Methods)

Section cross-reference(s): 10

ST **bacteria** enumeration **fluorescence** luminescence
luciferase GFP app

IT Luminescence spectroscopy
(bioluminescence; sensitivity of detection of **bacteria** with
fluorescent and luminescent phenotypes using different instruments)

IT Proteins, specific or class

RL: ANT (Analyte); ARU (Analytical role, unclassified); BPR (Biological
process); BSU (Biological study, unclassified); BUU (Biological use,
unclassified); **ANST (Analytical study)**; BIOL (Biological study);
PROC (Process); USES (Uses)

(green fluorescent; sensitivity of detection of **bacteria** with
fluorescent and luminescent phenotypes using different instruments)

IT Apparatus

(luminometer; sensitivity of detection of **bacteria** with
fluorescent and luminescent phenotypes using different instruments)

IT **Bacteria** (Eubacteria)

CCD cameras

Escherichia coli

Fluorescence

Fluorometers

Fluorometry

Luminescence

Luminescence, bioluminescence

Luminescence spectroscopy

Pseudomonas putida

Salmonella enteritidis

(sensitivity of detection of **bacteria** with fluorescent and
luminescent phenotypes using different instruments)

IT 9014-00-0, Bacterial luciferase 61970-00-1, Firefly Luciferase
RL: ANT (Analyte); ARU (Analytical role, unclassified); BPR (Biological
process); BSU (Biological study, unclassified); BUU (Biological use,
unclassified); **ANST (Analytical study)**; BIOL (Biological study);
PROC (Process); USES (Uses)

(sensitivity of detection of **bacteria** with fluorescent and
luminescent phenotypes using different instruments)

RE.CNT 9 THERE ARE 9 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 47 OF 179 CA COPYRIGHT 2003 ACS

AN 134:112423 CA

TI Multispectral bacterial identification

AU Tanner, Michael A.; Coleman, William J.; Everett, Christine L.; Robles,
Steven J.; Dilworth, Michael R.; Yang, Mary M.; Youvan, Douglas C.

CS Kairos Scientific, Inc., Santa Clara, CA, USA

SO Proceedings of SPIE-The International Society for Optical Engineering
(2000), 3913(In-Vitro Diagnostic Instrumentation), 45-53
CODEN: PSISDG; ISSN: 0277-786X

PB SPIE-The International Society for Optical Engineering

DT Journal

LA English

AB A multi spectral optical technique was developed to simultaneously
classify individual bacterial cells within mixed populations. Multi
spectral Bacterial Identification (mBID) combines innovations in
microscopy with a software anal. program to measure and characterize the
fluorescence signals from multiplexed 16S rRNA probes hybridized
to populations of different **bacteria**. Software was developed to
identify individual **bacteria** at the level of species within
these mixed populations. TO test the feasibility of mBID, we examd. the
fluorescence emissions from a mixt. of probes specific

for individual species of known **bacteria** from the American Type Culture Collection. Currently, up to seven species can be detected simultaneously by **fluorescence** microscopy. An eighth signal was reserved for a universal probe to control for **fluorescence** intensity. MBID can also be used to identify uncultured microorganisms. We plan to couple this new multi spectral technol. to existing identification technologies that utilize 16S rRNA sequence alignment. Using this integrated identification protocol, **bacteria** that may be assocd. with chronic conditions will be identified first by analyzing their 16S rDNA sequences and then by visualizing them with fluorescent probes hybridized to their 16S rRNA in situ.

- CC 9-4 (Biochemical Methods)
 Section cross-reference(s): 10, 14
- ST multispectral bacterial identification **fluorescence** microscopy
 rRNA specific probe vaginosis
- IT rRNA
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (16 S; multispectral bacterial identification)
- IT Vagina
 (disease; multispectral bacterial identification)
- IT Arthrobacter oxidans
 Bacillus subtilis
Bacteria (Eubacteria)
 Corynebacterium flavescens
 Escherichia coli
 Flexibacter maritimus
Fluorescence microscopy
 Lactobacillus delbrueckii lactis
 Leptothrix discophora
 Nucleic acid hybridization
 (multispectral bacterial identification)
- IT 312978-93-1D, conjugate with Bodipy 564/570
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (Arthrobacter globiformis 16S rRNA specific probe; multispectral bacterial identification)
- IT 312978-92-0D, conjugate with Bodipy R6G
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (Bacillus subtilis 16S rRNA specific probe; multispectral bacterial identification)
- IT 312978-95-3D, conjugate with Cy5
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (Corynebacterium flavescens 16S rRNA specific probe; multispectral bacterial identification)
- IT 312978-91-9D, conjugate with Pacific Blue
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (Escherichia coli 16S rRNA specific probe; multispectral bacterial identification)
- IT 312978-96-4D, conjugate with Cy5.5
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (Flexibacter maritimus 16S rRNA specific probe; multispectral bacterial identification)
- IT 312978-90-8D, conjugate with Alexa 350
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (Lactobacillus lactis 16S rRNA specific probe; multispectral bacterial identification)
- IT 312978-94-2D, conjugate with Bodipy 581/591
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES

(Uses)
 (Leptothrix discophora 16S rRNA specific probe; multispectral bacterial identification)

IT 320797-24-8D, conjugate with Bodipy 493/503
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (bacterial 16S rRNA specific probe; multispectral bacterial identification)

IT 989-38-8D, conjugate with oligonucleotide probe 121207-31-6D, Bodipy 493/503, conjugate with oligonucleotide probe 144377-05-9D, conjugate with oligonucleotide probe 150152-69-5D, Bodipy 581/591, conjugate with oligonucleotide probe 150173-89-0D, Bodipy 564/570, conjugate with oligonucleotide probe 172777-84-3D, Cy5.5, conjugate with oligonucleotide probe 215868-31-8D, Pacific Blue, conjugate with oligonucleotide probe 244636-14-4D, Alexa 350, conjugate with oligonucleotide probe
 RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
 (Uses)
 (multispectral bacterial identification)

RE.CNT 28 THERE ARE 28 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 50 OF 179 CA COPYRIGHT 2003 ACS
 AN 133:198214 CA
 TI Use of **fluorescence** for characterizing source and speciation of aquatic humic substances
 AU McKnight, Diane M.; Klapper, Lisa; Hood, Eran W.; Boyer, Elizabeth W.
 CS Institute of Arctic and Alpine Research, University of Colorado, Boulder, CO, 80309, USA
 SO Preprints of Extended Abstracts presented at the ACS National Meeting, American Chemical Society, Division of Environmental Chemistry (2000), 40(2), 659-660
 CODEN: PEACF2; ISSN: 1524-6434
 PB American Chemical Society, Division of Environmental Chemistry
 DT Journal
 LA English
 AB Aq. humic substances are defined as a heterogeneous class of moderate mol. wt., yellow-colored bio-mols. which are present in all natural water. Spectroscopic characterization of humic substances provides only limited structural information due to the heterogeneity within fulvic acids, although spectroscopic characterization is useful in quantifying differences among fulvic acids. **Fluorescence**, which can be detected at fulvic acid concns. present in most natural water, potentially contains information about the source and speciation of humic substances. Humic substances derived from degrdn. of **microbial** matter have distinctive **fluorescence** characteristics vs. those derived from plant and soil matter. Based on differences in representative **excitation emission** matrixes, a simple index was developed which could be used in a field study involving a large no. of samples. This index is the ratio of **emission** at 450 nm to the **emission** at 500 nm for an **excitation** of 370 nm. This index has a value of .apprx.1.8 for **microbially**-derived fulvic acids and .apprx.1.25 for terrestrially-derived fulvic acids. Since other chem. characteristics of fulvic acid vary between these 2 sources of org. matter, these results may indicate a significant seasonal change in fulvic acid reactivity. Thus, **fluorescence** measurements may provide a tool to est. humic substance reactivity in different environments and provide ancillary data to understand C cycling in aq. ecosystems.

CC 61-3 (Water)
 Section cross-reference(s): 80
 ST aquatic humic substance speciation source characterization;
fluorescence characterization aquatic humic substance
 IT Fulvic acids
 Humus

RL: ANT (Analyte); OCU (Occurrence, unclassified); **ANST (Analytical study)**; OCCU (Occurrence)
(aq.; using **fluorescence** to characterize sources and speciation of aquatic humic substances)

IT **Fluorescence**
(using **fluorescence** to characterize sources and speciation of aquatic humic substances)

IT 7732-18-5, Water, analysis

RL: AMX (Analytical matrix); **ANST (Analytical study)**
(natural water; using **fluorescence** to characterize sources and speciation of aquatic humic substances)

RE.CNT 9 THERE ARE 9 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 56 OF 179 CA COPYRIGHT 2003 ACS

AN 132:185144 CA

TI A novel method for detection of viable Giardia cysts in water samples

AU Jarmey-Swan, C.; Gibbs, R. A.; Ho, G. E.; Bailey, I. W.; Howgrave-Graham, A. R.

CS Analytical Services, Umgeni Water, Pietermaritzburg, 3200, S. Afr.

SO Water Research (2000), 34(6), 1948-1951

CODEN: WATRAG; ISSN: 0043-1354

PB Elsevier Science Ltd.

DT Journal

LA English

AB Assessing Giardia viability is a major requirement for public health purveyors and the water industry. Several indicators of viability (e.g., stains, excystation, animal infectivity) have been used to enumerate cysts with varying degrees of success. A combined detection-viability method for use in water samples would be useful to detect and det. cyst viability in source and drinking water and disinfection efficacy at treatment plants. Distd. water samples were seeded with purified Giardia cysts and incubated with fluorescein diacetate (FDA), initially to stain viable cysts, followed by tetra-Methyl red labeled anti-Giardia monoclonal antibodies (TMR) for confirmation of identity. As a result of FDA staining, green **fluorescence** of intact viable cysts was obsd. microscopically using a 450-490 nm exciter filter; non-viable cysts were not stained. Giardia cysts reacted pos. with TMR and glowed red using a triple band microscope filter with **excitations** of 400/450/570 nm. At this wavelength, a combination of FDA and TMR stained viable cysts green internally with a red wall while non-viable cysts only stained red. This simple, reliable, quick method allowed differentiation of Giardia cysts in water samples while simultaneously detg. their viability.

CC 61-3 (Water)

Section cross-reference(s): 10, 80

ST Giardia cyst viability detn water; staining monoclonal antibody identification viability confirmation Giardia cyst; fluorescein diacetate staining Giardia cyst water; tetramethylrhodamine labeled monoclonal antibody Giardia cyst identification

IT Cyst, **microbial**

(Giardia; staining/monoclonal antibody method to simultaneously detect, identify, and confirm viability of Giardia cysts in water)

IT Giardia

(cysts; staining/monoclonal antibody method to simultaneously detect, identify, and confirm viability of Giardia cysts in water)

IT Antibodies

RL: ARG (Analytical reagent use); NUU (Other use, unclassified); **ANST (Analytical study)**; USES (Uses)

(monoclonal, tetramethylrhodamine-labeled; staining/monoclonal antibody method to simultaneously detect, identify, and confirm viability of Giardia cysts in water)

IT 596-09-8, Fluorescein diacetate

RL: ARG (Analytical reagent use); NUU (Other use, unclassified); **ANST (Analytical study)**; USES (Uses)

(cyst staining with; staining/monoclonal antibody method to simultaneously detect, identify, and confirm viability of Giardia cysts in water)

IT 70281-37-7, Tetramethylrhodamine
 RL: ARG (Analytical reagent use); NUU (Other use, unclassified); ANST (Analytical study); USES (Uses)
 (monoclonal antibodies labeled with; staining/monoclonal antibody method to simultaneously detect, identify, and confirm viability of Giardia cysts in water)

IT 7732-18-5, Water, analysis
 RL: AMX (Analytical matrix); ANST (Analytical study)
 (source water and drinking water; staining/monoclonal antibody method to simultaneously detect, identify, and confirm viability of Giardia cysts in water)

RE.CNT 9 THERE ARE 9 CITED REFERENCES AVAILABLE FOR THIS RECORD
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 60 OF 179 CA COPYRIGHT 2003 ACS
 AN 131:141627 CA
 TI **Fluorescence** of dipicolinic acid as a possible component of the observed UV **emission** spectra of bacterial **spores**
 AU Nudelman, Raphael; Feay, Nicole; Hirsch, Mathew; Efrima, Shlomo; Bronk, Burt
 CS Mantech Environmental Technology Inc., USA
 SO Proceedings of SPIE-The International Society for Optical Engineering (1999), 3533(Air Monitoring and Detection of Chemical and Biological Agents), 190-195
 CODEN: PSISDG; ISSN: 0277-786X
 PB SPIE-The International Society for Optical Engineering
 DT Journal
 LA English
 AB Dipicolinic acid (DPA) and the Ca²⁺ complex of DPA (CaDPA) are well-known and are major chem. components of bacterial **spores**. DPA's native **fluorescence** is very weak and is thought to be completely masked by the **fluorescence** of tryptophan when this compd. is present. Thus **fluorescence** related to DPA in **spores** is assumed by many authors to be completely absent. We show that the **fluorescence** of CaDPA is substantial for **excitation** between about 290 nm and 310 nm with **emission** peaking near 400 nm. This **emission** is at the long wavelength tail for **emission** from tryptophan. We examine whether the **emission** of CaDPA could contribute to the total **emission** spectrum when bacterial **spores** are present in an aerosol, for **excitation** wavelengths in the neighborhood of 310 nm. In this report we present measurements of **fluorescence excitation** and **emission** for CaDPA and compare them with that of DPA and tryptophan.

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 4, 10
 ST bacterium **spore** dipicolinate fluorometry
 IT **Bacteria** (Eubacteria)
 Environmental analysis
 Fluorometry
Spore
 (fluorescence of dipicolinic acid as a possible component of obsd. UV **emission** spectra of bacterial **spores**)

IT 73-22-3, L-Tryptophan, analysis 499-83-2, Dipicolinic acid 499-83-2D, Dipicolinic acid, calcium complex 7440-70-2D, Calcium, complex with dipicolinic acid, analysis
 RL: ANT (Analyte); BOC (Biological occurrence); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); OCCU (Occurrence)
 (fluorescence of dipicolinic acid as a possible component of obsd. UV **emission** spectra of bacterial **spores**)

RE.CNT 15 THERE ARE 15 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 64 OF 179 CA COPYRIGHT 2003 ACS
AN 130:278941 CA
TI Fluorescent biological particle detection system
IN Ho, Jim Yew-wah
PA Her Majesty the Queen In Right of Canada, as Represented by the Minister
of, Can.
SO U.S., 18 pp., Cont.-in-part of U.S. 5,701,012.
CODEN: USXXAM
DT Patent
LA English
FAN.CNT 2

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	US 5895922	A	19990420	US 1997-863023	19970523
	US 5701012	A	19971223	US 1996-616475	19960319
PRAI	US 1996-616475		19960319		

AB A process and app. are provided for detection of viable and potentially hazardous biol. particles which may be dispersed in a particulate-contg. airstream. The process comprises directing each of the contained particles along a linear path through air, in a sequential manner, and sampling them for detn. of their size, whether they are biol. and viable, and whether they are present in concns. greater than background levels. The particle size identifies the particles as respirable or not and the particles are characterized as biol. and viable by subjecting each particle in turn, to 340 nm, UV laser light and looking for the **emission of fluorescence** which is typically emitted from **bacteria** or bacterial **spore**. **Fluorescence** detected in the 400-540 nm range signals the presence of NADH, which is indicative of biol. activity or viability. A more compact, and power-saving app. results with the preferential use of a solid state, UV laser, which is actuated only when the particle is passing the laser and only if it is deemed to be a biol. viable candidate.

IC ICM G01N021-64
NCL 250491200
CC 9-5 (Biochemical Methods)
ST **fluorescence** biol particle detection system
IT Particles
(Biol.; fluorescent biol. particle detection system)
IT Air analysis
Apparatus
Bacteria (Eubacteria)
Biochemical molecules
Fluorometry
Particle size
Sampling
Spore
UV lasers
(fluorescent biol. particle detection system)

IT 58-68-4, NADH
RL: ANT (Analyte); BSU (Biological study, unclassified); **ANST**
(**Analytical study**); BIOL (Biological study)
(fluorescent biol. particle detection system)

RE.CNT 10 THERE ARE 10 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 67 OF 179 CA COPYRIGHT 2003 ACS
AN 130:179461 CA
TI Two-photon **excitation** in **fluorescence** lifetime imaging
AU Gerritsen, Hans C.; Vroom, Jurrien; Sytsma, Joost
CS Debye Institute, Utrecht University, Utrecht, 3508 TA, Neth.
SO Fluorescence Microscopy and Fluorescent Probes, [Based on the Proceedings

of the Conference on Fluorescence Microscopy and Fluorescent Probes], 2nd, Prague, Apr. 9-12, 1997 (1998), Meeting Date 1997, 55-62. Editor(s): Slavik, Jan. Publisher: Plenum, New York, N. Y.
CODEN: 67BTAH

DT Conference
LA English
AB The authors employ time-gated lifetime imaging implemented in a two-photon **excitation** scanning microscope. The examples given demonstrate that **fluorescence** lifetime contrast is particularly suitable for in-depth imaging expts. Lifetime imaging enables the discrimination of multiple probes based on their differences in lifetime. The technol. is applied to the imaging of biofilms and human skin using Acridine-Orange.
CC 9-4 (Biochemical Methods)
Section cross-reference(s): 6, 13
ST microscopy two photon **excitation fluorescence** lifetime imaging
IT Staining, biological
(fluorescent; two-photon **excitation in fluorescence** lifetime imaging)
IT Tooth
(plaque; two-photon **excitation in fluorescence** lifetime imaging)
IT Imaging
(time-gated lifetime; two-photon **excitation in fluorescence** lifetime imaging)
IT Biofilm **bacteria**
Fluorescent indicators
Fluorescent substances
Skin
(two-photon **excitation in fluorescence** lifetime imaging)
IT Microscopy
(two-photon **excitation** scanning; two-photon **excitation in fluorescence** lifetime imaging)
IT Laser induced **fluorescence**
Photoexcitation
(two-photon; two-photon **excitation in fluorescence** lifetime imaging)
IT 65-61-2, Acridine Orange
RL: ARU (Analytical role, unclassified); BSU (Biological study, unclassified); PRP (Properties); **ANST (Analytical study)**; BIOL (Biological study)
(two-photon **excitation in fluorescence** lifetime imaging)
RE.CNT 18 THERE ARE 18 CITED REFERENCES AVAILABLE FOR THIS RECORD
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L17 ANSWER 79 OF 179 CA COPYRIGHT 2003 ACS
AN 128:45574 CA
TI Fluorescent biological particle detection system
IN Ho, Jim Yew-wah
PA Her Majesty the Queen In Right of Canada, as Represented by the Minister of, Can.
SO U.S., 17 pp.
CODEN: USXXAM
DT Patent
LA English
FAN.CNT 2

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	US 5701012	A	19971223	US 1996-616475	19960319
	US 5895922	A	19990420	US 1997-863023	19970523
PRAI	US 1996-616475		19960319		
AB	A process and app. are provided for detection of viable and potentially				

hazardous biol. particles which may be dispersed in a particulate-contg. airstream. The process comprises directing each of the contained particles along a linear path through air, in a sequential manner, and sampling them for detn. of their size, whether they are biol. and viable, and whether they are present in concns. greater than background levels. The particle size identifies the particles as respirable or not and the particles are characterized as biol. and viable by subjecting each particle in turn, to 340 nm, UV laser light and looking for the **emission of fluorescence** which is typically emitted from **bacteria** or bacterial **spore**. **Fluorescence** detected in the 400-540 nm range signals the presence of NAD hydrogen, which is indicative of biol. activity or viability.

IC ICM G01N021-64
NCL 250461200
CC 9-5 (Biochemical Methods)
Section cross-reference(s): 10
ST **fluorescence** biol particle detection system
IT **Spore**
(Bacterial; fluorescent biol. particle detection system)
IT Particles
(Biol.; fluorescent biol. particle detection system)
IT Air analysis
Apparatus
Bacteria (Eubacteria)
Biochemical molecules
Fluorometers
Fluorometry
UV lasers
(fluorescent biol. particle detection system)
IT 58-68-4, NADh
RL: ANT (Analyte); **ANST (Analytical study)**
(fluorescent biol. particle detection system)

L17 ANSWER 88 OF 179 CA COPYRIGHT 2003 ACS
AN 127:80177 CA
TI Noninvasive monitoring of the physiological state of **microbial** cultures
AU Lestan, D.; Perdih, A.
CS Centre Soil Environmental Sci., Dep. Agronomy, Biotechnical Faculty, Univ. Ljubljana, Ljubljana, Slovenia
SO Acta Chimica Slovenica (1997), 44(1), 1-15
CODEN: ACSLE7; ISSN: 1318-0207
PB Slovenian Chemical Society
DT Journal; General Review
LA English
AB A review with 41 refs. In attempts to improve the performance of bioprocess modeling and control it is becoming clear that alternative methods for accessing information from biol. systems, better suited to the nature of living systems, have to be developed. Three considerably different approaches have been proposed to access this information and are reviewed here. The optical approach relies on enzymes or metabolites which change their optical absorption or **fluorescence emission** as a function of specific or induced cellular alterations. ³¹P NMR can be used for the detn. of energy-rich P compds. Online monitoring of the physiol. state of the living matter in bioreactors uses knowledge-based recognition systems to assess variables that indicate the physiol. state.

CC 16-0 (Fermentation and Bioindustrial Chemistry)
Section cross-reference(s): 9, 10
ST review microorganism culture noninvasive monitoring; optical spectroscopy **microbial** culture monitoring review; phosphorus ³¹ NMR culture monitoring review; recognition system **microbial** culture monitoring review
IT Computer application

(expert systems; noninvasive monitoring of physiol. state of **microbial** cultures)

IT Bioreactors
Biotechnology
Fluorescence
Microorganism
Optical absorption
Spectroscopy
UV and visible spectroscopy
(noninvasive monitoring of physiol. state of **microbial** cultures)

IT Energy-rich phosphates
RL: ANT (Analyte); **ANST (Analytical study)**
(noninvasive monitoring of physiol. state of **microbial** cultures)

IT NMR spectroscopy
(phosphorus-31; noninvasive monitoring of physiol. state of **microbial** cultures)

L17 ANSWER 98 OF 179 CA COPYRIGHT 2003 ACS
AN 124:64699 CA
TI Aerosol-**fluorescence** spectrum analyzer: real-time measurement of **emission** spectra of airborne biological particles
AU Hill, Steven C.; Pinnick, Ronald G.; Nachman, Paul; Chen, Gang; Chang, Richard K.; Mayo, Michael W.; Fernandez, Gilbert L.
CS Army Research Laboratory, White Sands Missile Range, NM, 88002-5501, USA
SO Applied Optics (1995), 34(30), 7149-55
CODEN: APOPAI; ISSN: 0003-6935
PB Optical Society of America
DT Journal
LA English
AB We have assembled an aerosol-**fluorescence** spectrum analyzer (AFS), which can measure the **fluorescence** spectra and elastic scattering of airborne particles as they flow through a laser beam. The aerosols traverse a scattering cell where they are illuminated with intense (50 kW/cm²) light inside the cavity of an argon-ion laser operating at 488 nm. This AFS can obtain **fluorescence** spectra of individual dye-doped polystyrene microspheres as small as 0.5 .mu.m in diam. The spectra obtained from microspheres doped with pink and green-yellow dyes are clearly different. We have also detected the **fluorescence** spectra of airborne particles (although not single particles) made from various biol. materials, e.g., Bacillus subtilis **spores**, Bacillus anthracis **spores**, riboflavin, and tree leaves. The AFS may be useful in detecting and characterizing airborne **bacteria** and other airborne particles of biol. origin.

CC 59-1 (Air Pollution and Industrial Hygiene)
ST airborne biol aerosol particle detn; bioaerosol detn air laser **fluorescence**

IT Air analysis
Particles
(real-time measurement of **emission** spectra of airborne biol. particles)

IT Bacillus anthracis
Bacillus subtilis
(**spores**; real-time measurement of **emission** spectra of airborne biol. particles)

IT Aerosols
(airborne, biol., real-time measurement of **emission** spectra of airborne biol. particles)

IT Air pollution
(particulate, real-time measurement of **emission** spectra of airborne biol. particles)

IT 83-88-5, Riboflavin, analysis
RL: ANT (Analyte); **ANST (Analytical study)**

(real-time measurement of **emission** spectra of airborne biol.
particles)

L17 ANSWER 100 OF 179 CA COPYRIGHT 2003 ACS

AN 124:22668 CA

TI Green fluorescent protein as a new expression marker in mycobacteria

AU Kremer, Laurent; Baulard, Alain; Estaquier, Jerome; Poulain-Godefroy,
Odile; Loch, Camille

CS Laboratoire de Microbiologie Genetique et Moleculaire, INSERM, Lille
Cedex, F-59019, Fr.

SO Molecular Microbiology (1995), 17(5), 913-22

CODEN: MOMIEE; ISSN: 0950-382X

PB Blackwell

DT Journal

LA English

AB This study describes the use and the advantages of the green fluorescent protein (GFP) as a reporter mol. for mycobacteria. The gfp gene from *Aequorea victoria* was placed under the control of the HSP60 promoter in the shuttle vector pGFM-11. The gfp expression in the recombinant *Mycobacterium smegmatis* and BCG was readily detected on agar plates by the development of an intense green **fluorescence** upon irradiation with long-wave UV light. In mycobacteria containing a pGFM-11 derivative that lacks the hsp60 promoter, no **fluorescence** was observed. However, this plasmid was successfully used as a promoter-probe vector to identify RCG promoters. The **fluorescence emission** of GFP in mycobacteria harboring pGFM-11 and grown in liquid media could be quantified by spectrofluorimetry. This allowed for easy assessment of drug susceptibility. As GFP does not require the addition of substrates or cofactors, the green fluorescent bacilli could be directly observed within infected macrophages using **fluorescence** and laser confocal microscopy, or in tissue sections of infected mice. Finally, infected cells or free-living recombinant mycobacteria could also be analyzed by flow cytometry. The GFP thus appears to be a convenient reporter for mycobacteria, allowing tracing of recombinant mycobacteria, isolation of promoters with interesting properties, in vivo drug testing, and the development of new diagnostic tools.

CC 3-1 (Biochemical Genetics)

Section cross-reference(s): 9, 10

ST green fluorescent protein expression marker mycobacteria

IT Gene, **microbial**

RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
(Uses)

(gfp, reporter; green fluorescent protein as an expression marker in mycobacteria)

IT *Aequorea victoria*

Mycobacterium

Mycobacterium BCG

Mycobacterium smegmatis

(green fluorescent protein as an expression marker in mycobacteria)

IT Plasmid and Episome

(pGFM-11; green fluorescent protein as an expression marker in mycobacteria)

IT Spectrochemical analysis

(fluorometric, green fluorescent protein as an expression marker in mycobacteria)

IT Proteins, specific or class

RL: ARG (Analytical reagent use); **ANST (Analytical study)**; USES
(Uses)

(green fluorescent, green fluorescent protein as an expression marker in mycobacteria)

IT Genetic element

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES
(Uses)

(promoter, hsp60; green fluorescent protein as an expression marker in

mycobacteria)

- L17 ANSWER 105 OF 179 CA COPYRIGHT 2003 ACS
AN 123:51436 CA
TI Spectroscopic properties of tryptophan and **bacteria**
AU Tang, G. C.; Yang, Y. L.; Huang, Z. Z.; Hua, W.; Zhou, F.; Cosloy, S.;
Alfano, R. R.
CS City College, City University of New York, New York, NY, 10031, USA
SO Proceedings of SPIE-The International Society for Optical Engineering
(1995), 2387, 169-72
CODEN: PSISDG; ISSN: 0277-786X
DT Journal
LA English
AB **Fluorescence** spectra of tryptophan and **bacteria** were
measured at different concns. using a Mediscience CD-Scan unit. The
emission spectra of tryptophan were obtained using an
excitation wavelength at 280 nm. The **excitation** spectra
were obtained at the **emission** of 340 nm. The min. detectable
concn. of tryptophan was 10⁻⁸ M. The **emission** spectra for
bacteria were probed at 340 nm. The min. detectable no. of
bacteria in a beam of the **excitation** light was detd. to
be about 30. Assuming that the **emission** band at 340 nm of
bacteria comes from tryptophan, the no. of tryptophan per
bacterium was estd. to be 108. This approach to det. the no. is almost
consistent with that obtained using a wt. method.
CC 9-5 (Biochemical Methods)
Section cross-reference(s): 10
ST **fluorescence bacteria** tryptophan detection
IT **Bacteria**
Escherichia coli
Fluorescence
(spectroscopic properties of tryptophan and **bacteria**)
IT 54-12-6, DL-Tryptophan
RL: ANT (Analyte); PRP (Properties); ANST (Analytical study)
(spectroscopic properties of tryptophan and **bacteria**)
- L17 ANSWER 110 OF 179 CA COPYRIGHT 2003 ACS
AN 121:200084 CA
TI In-vivo **fluorescence** detection and imaging of
porphyrin-producing **bacteria** in the human skin and in the oral
cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma
AU Konig, Karsten; Schneckenburger, Herbert; Hemmer, Joerg; Tromberg, Bruce;
Steiner, Rudolf
CS Beckman Laser Institute, Irvine, CA, 92715, USA
SO Proceedings of SPIE-The International Society for Optical Engineering
(1994), 2135(Advances in Laser and Light Spectroscopy to Diagnose Cancer
and Other Diseases), 129-38
CODEN: PSISDG; ISSN: 0277-786X
DT Journal
LA English
AB Certain **bacteria** are able to synthesize metal-free fluorescent
porphyrins and can therefore be detected by sensitive autofluorescence
measurements in the red spectral region. The porphyrin-producing
bacterium *Propionibacterium acnes*, which is involved in the pathogenesis
of acne vulgaris, was localized in human skin. Spectrally-resolved
fluorescence images of **bacteria** distribution in the face
were obtained by a slow-scan CCD camera combined with a tunable liq.
crystal filter. The structured autofluorescence of dental caries and
dental plaque in the red is caused by oral **bacteria**, like
Bacteroides or *Actinomyces odontolyticus*. "Caries images" were created by
time-gated imaging in the ns-region after ultrashort laser
excitation. Time-gated measurements allow the suppression of
backscattered light and non-porphyrin autofluorescence. Biopsies of oral
squamous cell carcinoma exhibited red autofluorescence in necrotic regions

and high concns. of the porphyrin-producing bacterium *Pseudomonas aeruginosa*. These studies suggest that the temporal and spectral characteristics of bacterial autofluorescence can be used in the diagnosis and treatment of a variety of diseases.

CC 9-5 (Biochemical Methods)

Section cross-reference(s): 8, 10, 14

ST **fluorescence** imaging porphyrin producing bacterium; acne vulgaris skin **fluorescence** imaging; dental caries bacterium **fluorescence** imaging; oral squamous cell carcinoma **fluorescence** imaging

IT Mouth

Skin

(in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT Porphyrins

RL: ANT (Analyte); BOC (Biological occurrence); BSU (Biological study, unclassified); MFM (Metabolic formation); **ANST (Analytical study)**; BIOL (Biological study); FORM (Formation, nonpreparative); OCCU (Occurrence)

(in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT **Bacteria**

(porphyrin-producing; in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT Tooth

(disease, caries, in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT Imaging

(fluorescent, in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT Spectrochemical analysis

(fluorometric, in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT Mouth

(neoplasm, squamous cell carcinoma, in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

IT Acne

(vulgaris, in-vivo **fluorescence** detection and imaging of porphyrin-producing **bacteria** in human skin and in oral cavity for diagnosis of acne vulgaris, caries, and squamous cell carcinoma)

L17 ANSWER 114 OF 179 CA COPYRIGHT 2003 ACS

AN 120:265052 CA

TI Online, non-destructive biomass determination of bacterial biofilms by fluorometry

AU Angell, Peter; Arrage, Andrew A.; Mittelman, Marc W.; White, David C.

CS Cent. Environ. Biotechnol., Knoxville, TN, 37932, USA

SO Journal of Microbiological Methods (1993), 18(4), 317-27

CODEN: JMIMDQ; ISSN: 0167-7012

DT Journal

LA English

AB The lack of online methodol. for the detn. of **microbial** biomass and activity of attached **bacteria** has severely limited the study of biofilm physiol. This study showed that the fluorescent

emission of arom. amino acids in **microbial** biofilms can be used to det. the biomass formed on 316 stainless steel coupons. Cells resuspended from the substratum were enumerated by viable and acridine orange counts showing correlation coeffs. of 0.77 and 0.98, resp., when compared to the tryptophan **fluorescence**. Substrata treated with a fluorescent epoxy coating (F-150) showed no **fluorescence** that could be attributed to the microorganisms. Bioluminescent **emission** of an actively growing bioluminescent bacterium, *Vibrio harveyi*, was correlated with acridine orange counts ($r^2 = 0.95$) and **fluorescence** ($r^2 = 0.93$). The results of these studies suggest that **fluorescence** measurements can be used to monitor **microbial** biomass assocd. with various substrata. Coupled with bioluminescence measurements, this method provides information on both biomass constituents and metabolic activity, and therefore possibly an indicator of sub-lethal toxicity.

CC 9-5 (Biochemical Methods)

Section cross-reference(s): 16, 60

ST **bacteria** biofilm biomass detn fluorometry

IT Films

(bacterial, biomass of, detn. of, by fluorometry)

IT **Bacteria**

(biofilms, biomass detn. in, by fluorometry)

IT Biomass

(detn. of, of **bacteria** biofilms by fluorometry)

IT 53-57-6, NADPH

RL: **ANST (Analytical study)**

(detn. of in **bacteria** biofilm biomass detn.)

IT 58-68-4, NADH 73-22-3, Tryptophan, analysis

RL: **ANT (Analyte); ANST (Analytical study)**

(detn. of, in **bacteria** biofilm biomass detn.)

L17 ANSWER 117 OF 179 CA COPYRIGHT 2003 ACS

AN 119:45156 CA

TI Determination of fluorescent substances in microorganism, and its use for determination of microorganism viability

IN Hirotsuji, Junji; Yoshimura, Yumiko; Sugimoto, Masuo; Nakatsugawa, Naoki; Oota, Naomi

PA Mitsubishi Electric Corp, Japan

SO Jpn. Kokai Tokkyo Koho, 20 pp.

CODEN: JKXXAF

DT Patent

LA Japanese

FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	JP 05111394	A2	19930507	JP 1991-277705	19911024
	JP 2947305	B2	19990913		
PRAI	JP 1991-277705		19911024		

AB The fluorescent substances such as coenzyme F420 (I) in microorganism is detd. by direct irradiation of the microorganism with an **excitation** wavelength. Based on the **fluorescence** data, the content and concn. of the fluorescent substances can be calcd. and the viability of the microorganism detd. The method is easy and fast, and does not need to disrupt the cells. Fluorescent detn. of I in methane **bacteria** and the viability of the **bacteria** was shown. The result was comparable to that with the prior art.

IC ICM C12Q001-02

CC 10-6 (Microbial, Algal, and Fungal Biochemistry)

ST **fluorescence** detn microorganism viability; substance

fluorescence detn microorganism

IT Fluorescent substances

(fluorescent detn. of, in microorganism for detn. of viability)

IT Microorganism

(fluorescent substances in, detn. of, for detn. of microorganism)

viability)
 IT **Bacteria**
 (methanogenic, F420 in, detn. of, for detn. of **bacteria**
 viability)
 IT 64885-97-8, Coenzyme F420
 RL: ANT (Analyte); **ANST (Analytical study)**
 (**fluorescence** detn. of, in methane **bacteria**)

L17 ANSWER 128 OF 179 CA COPYRIGHT 2003 ACS
 AN 115:178859 CA
 TI Method and reagents for detecting microorganisms
 IN Monget, Daniel
 PA Biomerieux S. A., Fr.
 SO Eur. Pat. Appl., 15 pp.
 CODEN: EPXXDW
 DT Patent
 LA French
 FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	EP 424293	A1	19910424	EP 1990-420453	19901018
	EP 424293	B1	19950412		
	R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE				
	FR 2653447	A1	19910426	FR 1989-14087	19891020
	FR 2653447	B1	19911227		
	AT 121138	E	19950415	AT 1990-420453	19901018
	ES 2071067	T3	19950616	ES 1990-420453	19901018
	CA 2028059	C	20010213	CA 1990-2028059	19901019
	US 5336600	A	19940809	US 1992-961625	19921016
PRAI	FR 1989-14087	A	19891020		
	US 1990-600919	B3	19901022		

OS MARPAT 115:178859
 AB Microorganisms are detected using an aq. reaction media contg. a C source, a N source, and a marker whose luminous **emission** is modified as a consequence of development of the microorganism in the reaction media. The marker is, e.g. I (R1, R4, R5 = H, F, Cl, Br, alkyl, alkoxy, CO2H, amide, cyano; R2, R3 = H, F, Cl, Br, alkyl, alkoxy, carboxylate, CO2H, amide, cyano; or R3R2 form an unsatd. ring; X - OH, amine) or its anionic form. Six different antibiotics at 2 different concns. were added to tubes contg. Mueller Hinton media, glucose, resorufin, and water (pH 7.3). *Escherichia coli*, *Pseudomonas aeruginosa*, and *Staphylococcus aureus* (10 cells/mL) were added and fluorescence was measured after 18 h incubation at 35.degree.. If the microorganism was sensitive to the antibiotic, **fluorescence** was maximal; resistant microorganisms gave total extinction of the fluorescence.

IC ICM C12Q001-04
 ICS C12Q001-18; C07D265-38
 CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10
 ST microorganism detection resorufin deriv **fluorescence**; antibiotic microorganism sensitivity **fluorescence** resorufin

IT *Candida albicans*
Citrobacter freundii
Escherichia coli
Pseudomonas aeruginosa
Staphylococcus aureus
Streptococcus faecalis
 (detection of, by **fluorescence** assay, resorufin as marker in)

IT Microorganism
 (detection of, by luminescence assay, resorufin or deriv. as marker in)

IT Blood analysis
 Body fluid
 Cerebrospinal fluid
 Cosmetics

Food analysis
 Pharmaceutical analysis
 Urine analysis
 (microorganism detection in, by luminescence assay, resorufin or deriv.
 as marker for)

IT Antibiotics
 Fungicides and Fungistats
 (microorganism sensitivity to, resorufin or deriv. in luminescence
 assay of)

IT Aeromonas hydrophila
 Candida tropicalis
 Klebsiella pneumoniae
 Proteus vulgaris
 Torulopsis glabrata
 Vibrio alginolyticus
 (response of, to different sugars, in **fluorescence** assay with
 resorufin as marker, microorganism identification in relation to)

IT Spectrochemical analysis
 (fluorometric, microorganism detection by, resorufin in)

IT **Bacteria**
 (intestinal, detection of, by **fluorescence** assay, resorufin
 as marker in)

IT Spectrochemical analysis
 (luminescence, microorganism detection by, resorufin or deriv. as
 marker in)

IT 635-78-9, Resorufin 635-78-9D, Resorufin, anions
 RL: **ANST (Analytical study)**
 (as luminescent marker for microorganism detection)

IT 101490-85-1
 RL: **ANST (Analytical study)**
 (as substrate for .beta.-glucuronidase of Escherichia Coli detection)

IT 50-99-7, Glucose, analysis 57-50-1, Saccharose, analysis 63-42-3,
 Lactose 69-79-4 99-20-7, Trehalose 528-50-7, Cellobiose 585-99-9,
 Melibiose
 RL: **ANST (Analytical study)**
 (microorganism response to, in microorganism identification by
 resorufin **fluorescence** assay)

IT 60-54-8, Tetracycline 61-33-6, properties 114-07-8, Erythromycin
 35607-66-0 61477-96-1
 RL: **ANST (Analytical study)**
 (microorganism sensitivity to, resorufin in **fluorescence**
 assay of)

IT 9001-45-0
 RL: **ANST (Analytical study)**
 (resorufin-glucuronide as substrate for, for Escherichia coli
 detection)

L17 ANSWER 141 OF 179 CA COPYRIGHT 2003 ACS
 AN 110:150195 CA
 TI **Fluorescence** photometric determination of the coenzyme F420 to
 monitor anaerobic effluent purification
 AU Kaiser, G.; Frenzel, S.; Mauch, W.
 CS Fachgeb. Zuckertechnol., Tech. Univ. Berlin, Berlin, D-1000/65, Fed. Rep.
 Ger.
 SO Zuckerindustrie (Berlin, Germany) (1988), 113(10), 868-72
 CODEN: ZUCKDI; ISSN: 0344-8657
 DT Journal
 LA German
 AB A method is proposed for the detn. of coenzyme F420, which is
 characteristic for CH₄-forming **bacteria**, making it possible to
 monitor anaerobic effluent purifn. by relatively simple means. The only
 equipment required is a slightly modified photometer or nephelometer, a
 table centrifuge, and a membrane filter. By detg. the
fluorescence intensity at 470 nm over a wide spectrum of

excitation wavelengths (340-450 nm), the optimum reaction conditions (e.g. pH, solvent, redox state) for fluorimetric detection have been identified and the specificity of the method for coenzyme F420 established. Since for normal plant control it is not necessary to measure the molar concn. of coenzyme F420, the anal. is reduced to the detn. of the relative **fluorescence**. The latter reflects the state of effluent fermn. Disturbances in the course of fermn. are quickly manifested by a decrease in **fluorescence**, allowing prompt correction of the effluent purifn. process.

CC 7-1 (Enzymes)
 Section cross-reference(s): 16
 ST coenzyme F420 detn **fluorescence** methanogen fermn
 IT Fermentation
 (with methanogens, monitoring of, **fluorescence** method for)
 IT **Bacteria**
 (methanogenic, fermn. with, monitoring of, **fluorescence** method for)
 IT 64885-97-8, Coenzyme F420
 RL: ANT (Analyte); **ANST (Analytical study)**
 (detn. of, by **fluorescence** method, for monitoring fermn. by methanogenic **bacteria**)

L17 ANSWER 149 OF 179 CA COPYRIGHT 2003 ACS
 AN 105:130360 CA
 TI Determination of a concentration of **bacteria** in a suspension
 IN Kosarev, N. V.; Puckkov, E. O.
 PA All-Union Scientific-Research Institute of Applied Microbiology, USSR
 SO U.S.S.R.
 From: Otkrytiya, Izobret. 1986, (18), 127.
 CODEN: URXXAF

DT Patent
 LA Russian

FAN.CNT 1

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	SU 1231077	A1	19860515	SU 1984-3719309	19840328
PRAI	SU 1984-3719309		19840328		

AB The concn. of **bacteria** in a suspension is detd. by prepg. a bacterial suspension and recording the optical parameters. The accuracy and sensitivity of anal. are increased by adding up to 0.03-0.05% Triton X100, .ltoreq.0.3-3.0 mmole tris-ethylenediaminetetraacetate (sic) at pH 7.2-7.5, and .ltoreq.4.0-7.0 .mu.mole ethidium bromide to the suspension and increase of the **fluorescence** intensity is measured in sample with **bacteria** as compared to the samples without **bacteria** at 540-620 nm with **excitation** at 260-320 or 420-530 nm, and the concn. of **bacteria** is calcd. from the increase of the **fluorescence** intensity.

IC ICM C12N001-00
 ICS C12Q001-00
 ICI C12Q001-00, C12R001-01

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10
 ST **bacteria** detn suspension fluorometry
 IT **Bacteria**
 Microorganism

(detn. of, in suspensions, by fluorometry)
 IT Spectrochemical analysis
 (fluorometric, for **bacteria** in suspensions)
 IT 1239-45-8 9002-93-1
 RL: **ANST (Analytical study)**
 (in **bacteria** detn. in suspension by fluorometry)

L17 ANSWER 155 OF 179 CA COPYRIGHT 2003 ACS
 AN 103:3115 CA

TI Identification of bacterial pathogens by laser excited
fluorescence
 AU Coburn, J. T.; Lytle, F. E.; Huber, D. M.
 CS Dep. Chem., Purdue Univ., West Lafayette, IN, 47907, USA
 SO Analytical Chemistry (1985), 57(8), 1669-73
 CODEN: ANCHAM; ISSN: 0003-2700
 DT Journal
 LA English
 AB The title rapid method relies on the extent of aminopeptidase hydrolysis of a series of nonfluorescent L-amino acid .beta.-naphthylamides to produce the highly fluorescent .beta.-naphthylamine by the pathogen of interest. The luminescence background was composed of Raman and Rayleigh scatter, fluorescent impurities in the buffer, .beta.-naphthylamine **fluorescence** due to substrate decomn., and **emission** of the biol. matrix. The blank levels were systematically examd. and reduced to levels which allow the measurement of the fluorophore in the 0.1-nM range. Thus, unambiguous identification of pathogens at the 50,000 cell/mL level was achieved. This corresponds to 2-3 orders of magnitude fewer cells than needed by other techniques. Identification of pathogens at this level will drastically reduce the cell growth period from 48 to 6 h and afford a more rapid turnaround time for bacterial identification.

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10
 ST **bacteria** pathogen identification laser fluorometry;
 aminopeptidase profiling **bacteria** identification
 IT Pseudomonas syringae phaseolicola
 Xanthomonas phaseoli
 (identification of, by aminopeptidase profiling with laser fluorometry)
 IT Spectrochemical analysis
 (fluorometric, laser-induced, for aminopeptidase substrates, in
bacteria pathogen profiling)
 IT Fluorometers
 (laser, time-resoln., for **bacteria** pathogen identification)
 IT **Bacteria**
 (pathogenic, identification of, by aminopeptidase profiling with laser
 fluorometry)

IT 716-94-9 720-82-1 729-24-8 732-84-3 732-85-4 740-57-8
 1259-69-4 3326-63-4 3326-64-5 4357-95-3 4420-88-6 7182-70-9
 7424-15-9 7424-16-0 14525-44-1 16037-15-3
 RL: RCT (Reactant); RACT (Reactant or reagent)
 (hydrolysis of, by aminopeptidase in **bacteria** pathogen
 identification by laser fluorometry)
 IT 9031-94-1
 RL: ANST (Analytical study)
 (in **bacteria** pathogen identification by laser fluorometry)

L17 ANSWER 156 OF 179 CA COPYRIGHT 2003 ACS
 AN 102:75146 CA
 TI A comprehensive method for the measurement of **fluorescence**
 lifetime in picoseconds with a phase fluorometer and its application for
 the determination of **excitation** energy transfer rate in
 light-harvesting pigment antenna of green **bacteria**
 AU Fetisova, Z. G.; Kharchenko, S. G.; Blagoveshchenskii, Yu. N.; Borisov, A.
 Yu.
 CS USSR
 SO Vestnik Moskovskogo Universiteta, Seriya 16: Biologiya (1984), (4), 56-9
 CODEN: VMUBDF; ISSN: 0137-0952
 DT Journal
 LA Russian
 AB A universal phase fluorometric method for measuring **fluorescence**
 lifetime in picoseconds is based on mixing 2 **fluorescence** beams.
 The photocathode of the phase fluorometer registers the superimposed 2
 light components: the desired short-lived and long-lived components of the
 same object in a suitable solvent. This method was specially modified for

measuring the **fluorescence** of photosynthetic pigments (in vivo). The 2-component systems were obtained from 2 solns. of Na fluorescein in 0.001M NaOH, placed in the adjacent cuvettes, 1 of which contained KI as a quencher. Math. equations are given to do the necessary calcs. The method was used to measure the rate of **excitation** energy transfer in light-harvesting pigment antenna of Chlorobium limicola with a photoreaction center P840. About 90% bacteriochlorophyll c with light-harvesting property participates in transfer of **excitation** energy transfer to bacteriochlorophyll a and accomplishes it in 20-60 ps with an efficiency of >95% under unsatisfactory light conditions (for photosynthesis).

CC 9-5 (Biochemical Methods)
 Section cross-reference(s): 10, 11
 ST **fluorescence** lifetime picosecond detn; energy transfer
 Chlorobium **fluorescence** picosecond
 IT Energy transfer
 (fluorescence, in photosynthetic systems of green
 bacteria, method for measurement of)
 IT Pigments, **microbial**
 (light-harvesting, **excitation** energy transfer rate detn. in)
 IT Chlorobium limicola
 (photosynthetic systems of, **excitation** energy transfer in,
 method for detn. of)
 IT Fluorometry
 (picosecond, of **excitation** energy transfer in
 light-harvesting pigment antenna of green bacteria)
 IT Photosynthetic systems
 (reaction center, P840, of Chlorobium limicola, **excitation**
 energy transfer in, method for detn. of)
 IT 17499-98-8 53986-51-9
 RL: **ANST (Analytical study)**
 (of Chlorobium limicola, method for **excitation** energy
 transfer rate detn. in relation to)

L17 ANSWER 169 OF 179 CA COPYRIGHT 2003 ACS

AN 93:40914 CA

TI Two-parameter analysis of **microbial** cell constituents

AU Hutter, K. J.; Stoehr, M.

CS Inst. Exp. Pathol., Dtsch. Krebsforschungszent., Heidelberg, 6900, Fed.
 Rep. Ger.

SO Microbios Letters (1979), 10(39-40), 121-8

CODEN: MILEDM; ISSN: 0307-5494

DT Journal

LA English

AB Flow cytometry is a new assay to investigate different cellular constituents, e.g., the DNA and protein content of a wide variety of biol. specimens. Investigations have involved quant. fluorescent staining of monodisperse cell populations in liq. suspension at a flow rate of 1000 cells/s. The development of a new dual laser beam **excitation** device for flow cytometry anal. has enabled simultaneous measurement of the DNA and protein content of baker's yeast cells. The combination of light sources consisted of an Ar ion laser with an **excitation** wavelength of 488 nm, and a 2nd Ar ion laser at 360 nm wavelength. The **microbial** DNA content was stained by 4',6-diamidino-2-phenylindol while the protein content was fluorochromized by Sulforhodamine 101. This staining technique avoided any RNA digestion with pepsin in order to eliminate non-specific cytoplasm **fluorescence** and revealed a no. of cells in various phases of their life cycle. A correlation between the replicative and metabolic activity was visualized.

CC 9-6 (Biochemical Methods)

ST protein yeast flow cytometry; DNA yeast flow cytometry

IT Saccharomyces cerevisiae

(DNA and protein detn. in, by flow cytometry)

IT Deoxyribonucleic acids

Proteins

RL: ANT (Analyte); ANST (Analytical study)
(detn. of, in yeast by flow cytometry)

IT 28718-90-3

RL: ANST (Analytical study)
(DNA in yeast staining with, for detn. by flow cytometry)

IT 60311-02-6

RL: ANST (Analytical study)
(protein staining in yeast with, for detn. by flow cytometry)

L17 ANSWER 175 OF 179 CA COPYRIGHT 2003 ACS

AN 78:56378 CA

TI Method for simultaneous determination of histidine and histamine in biological liquids. Application to wines

AU Plumas, B.; Sautier, C.

CS Cent. Rech. Diet., Hop. Bichat, Fr.

SO Annales des Falsifications et de l'Expertise Chimique (1972), 65(703), 322-36

CODEN: AFECAT; ISSN: 0003-4274

DT Journal

LA French

AB Histidine (I) is sepd. from histamine (II) on an Amberlite CG-50 resin; II is eluted from the resin with 3N HCl. To det. I, the pH of the effluent is adjusted to 12.15 with NaOH, o-phthalic dialdehyde (C₆H₄(CHO)₂) is added, and after exactly 4 min the condensation reaction is stopped by the addn. of 3N HCl. A **fluorescence** that is stable for 30 min is obtained; activation 360 m.mu., and **emission** 450 m.mu.. To det. II, the pH is adjusted to 12.45. The precision and sensitivity of the method for II is 0.05 mg/l. and 0.025 .mu.g/ml, resp., and for I is 0.7 mg/l. and 0.75 .mu.g/ml. Anal. of about 50 wines of different origins gave higher values for I in white wine, 19.43-31.28 mg/l., than for red wines, 16.99-17.22 mg/l. White wines had less II, 2.69-4.90 mg/l., than red wines, 6.36-6.49 mg/l. In the course of malolactic reversion, the amt. of I decreased and II increased. This may be due to a simultaneous decarboxylation of I and malic acid. To avoid formation of II, **bacteria** should be selected that can decarboxylate malic acid without acting on I. 9 refs.

CC 16-1 (Fermentations)

ST wine histidine histamine detn; malate wine histamine; fluorometry histamine

IT Wine analysis

(histamine and histidine simultaneous detn. in)

IT 71-00-1, analysis

RL: ANT (Analyte); ANST (Analytical study)
(detn. of, in histamine presence)

IT 51-45-6, analysis

RL: ANT (Analyte); ANST (Analytical study)
(detn. of, in histidine presence)

L17 ANSWER 179 OF 179 CA COPYRIGHT 2003 ACS

AN 72:107693 CA

TI Direct fluorometric determination of bacterial nucleic acids

AU Launay, Bernard; Truhaut, Rene

CS Centre Rech. Toxicol., Fac. Pharm., Paris, Fr.

SO Comptes Rendus des Seances de l'Academie des Sciences, Serie D: Sciences Naturelles (1969), 269(25), 2614-17

CODEN: CHDDAT; ISSN: 0567-655X

DT Journal

LA French

AB Cultures of *Proteus vulgaris* and *Escherichia coli* were collected, washed with NaCl 0.16 + citrate 0.01M, and .apprx.20 g wet wt. of each bacterium was suspended in 20 ml of the washing soln. The **bacteria** were lysed with 0.25-0.5 ml of 20% Na dodecylsulfate in EtOH, and the lysate was diluted 200-fold with Tris-HCl buffer at pH 7.5. **Fluorescence**

was detd. in the presence of 8 $\mu\text{g/ml}$ of ethidium bromohydrate with **excitation** at 3400 \AA and **emission** at 5950 \AA . The increase in the **fluorescence** of ethidium bromohydrate was proportional to the vol. of the lysate. The method allowed the detn. of DNA at concns. $\geq 10^{-3} \mu\text{g/ml}$, corresponding to a level of $\approx 5 \times 10^4$ **bacteria/ml**.

CC 6 (Biochemical Methods)
ST fluorometry nucleic acids; nucleic acids fluorometry
IT Escherichia coli
(deoxyribonucleic acid detn. in)
IT Nucleic acids, deoxyribo-
RL: ANT (Analyte); **ANST (Analytical study)**
(detn. of, in **bacteria**)
IT Proteus
(vulgaris, deoxyribonucleic acid detn. in)

=>